

Blood Bank Refrigerator (Double door)



Introduction:

Blood bank refrigerator series products are applied to places needing storage at constant temperature, such as hospitals, epidemic prevention stations, blood banks, scientific research institutions, colleges and universities, biological pharmacies, genetic engineering, etc. The storing temperature in the body can be regulated by the temperature button on the control panel. The application is convenient and the performance is reliable.

Advantage:

- 1. Multiple sensors design: make sure the precise monitoring and programming.
- 2. Electric heated glass semi-foamed door, easy to observe, also keep a better cooling temperature.
- 3. Forced air refrigeration system, ensure inner temperature constant and uniform.
- 4. Adjustable shelves, improve space utilization; blood storage baskets can be classified.
- 5. Individual transparent inner door, effectively reduce cool air loss when opening the door.

Technical Parameters

Model		BBR-4V626	BBR-4V966
Capacity		626L	966L
Temp. Range		4±2℃	
Temp. Accuracy		0.1℃	
Control System		Microprocessor control	
Display		LED display	
Alarm		Audio and visual alarm for: High or low temperature, Sensor & system failure; Door ajar, Power failure, Low battery	
Refrigeration Type		Forced air refrigeration system	
Refrigerant		R134a, CFC-Free	
Defrost		Frostless Design	
Construction		Unibody Design, Mono-assembly foaming (Rigid polyurethane insulation material)	
Material	Internal Material	Stainless steel	
	External Material	Cold-rolled steel with anti-bacteria powder coating	
Door		Electric heated glass semi-foamed door with lock, Heating film preventing frosting, Sealing joint strip with magnetic force. Transparent plastic inner door*6	
Shelves		6*2, Adjustable shelves	
Casters		4 casters	
Consumption		529W	529W
Power Supply		AC110V/220V±10%, 50/60Hz	
Blood Bags (450ml)		10*36=360 pcs	16*36=576 pcs
Standard Accessory		Temperature printer, 12 shelves, 36 baskets, LED lamp, probe access port, door key	
Optional Accessory		Temperature recorder	
External Size (W*D*H) mm		1220*665*1885	1220*895*1885
Package Size (W*D*H) mm		1345*770*2020	1345*975*2020
Gross Weight		255kg	329kg