

Table Top Autoclave Class B Series

Autoclaves are equipment for sterilizing widely used in microbiology, medicine, podiatry, tattooing, body piercing, veterinary science, mycology, dentistry, and prosthetics fabrication. They vary in size and function depending on the media to be sterilized.



Features

- 1.Built-in open type water tank
- 2.Brand-new operation interface
- The LCD screen can display temperature, pressure, time, operating status, error and alarm etc.
- 3.Multiple security protective device
- Over temperature automatic protective device, multiple control and protection for steam generator, door safety interlock protection and dual chamber overpressure protection.
- 4.Efficient vacuum
- Efficient vacuum system, low working noise, vacuum lowest limit can reach over -90KPa.

Technical Parameter

Model	BKM-Z18A	BKM-Z23B
Capacity	18L	23L
Sterilization Class	Class B	
Sterilization Temp.	121℃/134℃	
Special Program	HIV, HBV long-sterilization program	
Drying System	Vacuum drying system	
Display	LCD display	
Control Precision	Temperature: 1℃ Pressure: 0.1bar	
Sterilization Data Output	Optional Mini printer Standard USB port	
Safety System	Pressure lock system Relief valve in case of over pressure Pressure or temperature over load protection Alarm for system failure, finish reminding, water level warning	
Water Supply System	Build-in water tank, easy to clean	
Water Tank Capacity	3.5L	
Tray Holder	3 pcs SS trays on SS shelf	
Chamber	Medical stainless steel 304 unibody structure	
Ambient Temp.	5~40℃	
Noise	≤65dB	
Power Supply	110/220V±10%, 60/50Hz	
Power Consumption	2000W	2000W
Chamber Size(mm)	φ247*360	φ247*470
External Size(mm)	600*495*410	700*495*410
Packing Size(mm)	810*610*590	810*610*590
Gross Weight(kg)	63	65